





Rocket® Seldinger Chest Drainage Sets




The Rocket® Seldinger Chest Drainage Set is indicated for the percutaneous insertion of a pleural drain for fluid and air and is particularly suitable for: ¹

- Pneumothorax, spontaneous or iatrogenic
- Effusions, benign or malignant, except those with high protein content
- Low viscosity empyema

DESCRIPTION	SIZE	ORDER CODE
Rocket® Seldinger Chest Drainage Kit with Safety Guard and fixation set Skin knife, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir-tree & wide-bore luer connectors, skin fixation set. Sterile, for single use in cartons of 5 units.	12FG	R54544-12-SG
	18FG	R54544-18-SG
Rocket® Seldinger Chest Drainage Tray Complete Ward Procedure Set Skin knife, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir -tree & wide-bore luer connectors, skin fi xation set, ward procedure pack including drapes, swabs, sponges, 10ml syringes, gallipot. Sterile, for single use in cartons of 5 units.	12FG	R54544-12-PK
	16FG	R54549-16-PK
	18FG	R54544-18-PK
	20FG	R54549-20-PK

TECHNICAL INFORMATION

- Improved patient comfort using small bore catheterisation^{2,3}
- Reliable insertion via Seldinger approach under direct vision, x-ray screening or ultrasound guidance
- 10cm dilators with safety stops to minimise the risks of over-insertion associated with the traditional 20cm dilator.



Drapes and a gown are included in all PK coded packs in response to NPSA/2008/RRR3

References:

- Clinical Evaluation of Seldinger Drainage Kit. Beaumont A C et al Proceedings of Intensive Care Society Meeting. December 1999.
- 2. Catheters smaller than 24 French gauge can be used for chest drains. Taylor PM. BMJ (1997) 315:186
- 3. Management of Malignant Pleural Effusions Antunes G, Neville E,Thorax (2000) 55:981-983